

Mio Medical Ltd

Jl Kyai Maja 43 Jakarta 12120

phone /fax : +6282113355756

sales@miomedical.com

Patient Monitor



New 2009 manufactured, GE Dash 5000 Patient Monitor Price US\$5,500

Features

- 12.1" Color Display
- NIBP - Dinamap Superstat
- IBP - 4 Channels
- 12 Lead ECG - 12SL, ACI-TIPI, Intellirate, AVOA Plus
- Spo2 - Masimo SET
- CO2 - Includes Capnoflex for Sidestream CO2 monitoring
- LED Alarm Handle
- Printer/Recorder
- 2 Batteries
- Wi-Fi®: IEEE 802.11 b/g
- New accessories as follows
- 5 or 12 Lead ECG Cable
- NIBP Cuff and Hose
- Spo2 Extension Cable and Finger Sensor
- CO2 Module
- Size Dash 5000 – 12.1 in.
- Type Active-matrix color LCD
- Resolution: Dash 3000 and 4000: 640 by 480 dpi, Dash 5000: 800 by 600 dpi
- Number of traces 7 (maximum)
- Number of seconds/trace 4.9 at 25 mm/sec
- Sweep speed 6.25, 12.5, 25 mm/sec (with erase bar)
- Controls
- TrimKnob control
- Five hard keys Standard Silence Alarm, Print, NBP Go/Stop, Zero All and Power On/Off. Dash 5000 adds Trends, NBP Auto, Admit/Discharge, Standby and Main View

- Remote control option Available
- Alarms
- Categories Patient status and system status
- Priority 4 levels – Crisis, Warning, Advisory and Message
- Notification Audible and visual
- Setting Default and individual
- Silencing 1 minute, current alarm only
- Pause 5 minutes in Adult ICU mode, 3 minutes in Neonatal ICU mode and 5 minute, 15 minute, or permanent pause in OR mode
- Volume Default 70 dB measured at 1 meter
- Invasive blood pressure
- Number of channels 1 to 4 (optional)
- Transducer sites Arterial, femoral artery, pulmonary arterial, central venous, right atrial, left atrial, intracranial and special
- Transducer requirements Excitation voltage: 5 V dc \pm 0.1%
- Transducer output 5 μ V/V/mmHg
- Input specifications
- Range -25 mmHg to 300 mmHg
- Offset \pm 150 mmHg
- Output specifications
- Frequency response dc to 50 Hz (-3 dB)
- Zero balance range \pm 150 mmHg
- Zero balance accuracy \pm 1 mmHg
- Zero balance drift \pm 1 mmHg over 24 hours
- Accuracy \pm 2% or \pm 1 mmHg, whichever is greater (exclusive of transducer)
- Alarms User-selectable upper and lower limits for systolic, diastolic and mean pressures

New 2009 manufactured, GE Dash 4000 Patient Monitor Price US\$2,500



- One internal battery
- NIBP Blood pressure adult cuff with hose (\$120 Value)
- SPO2 Pulse Ox preamp cable with Adult Finger Sensor (\$180 Value)
- ECG cable with leads (\$100 Value)
- Hospital Grade Power Cable
- PDF Operators Manual
- Noteworthy GE Dash 4000 Monitor Features COLOR LCD display
- ECG, SPO2, NIBP, RESP, and TEMP parameters
- TWO INVASIVE PRESSURE CHANNELS
- Built in Printer/Recorder

- Battery Operated
- 100-240volt, 50-60hz power flexibility



PhysioControl Lifepak 15 Patient Monitor
Defibrillator
Compare price: \$37,999.00
Our price: \$12900



GE Dash 5000 Patient Monitor
Compare price: \$12,999.00
Our price: \$7,499.00



GE Datex Ohmeda Cardiocap/5 Anesthesia
Monitor
Compare price: \$9,999.00
Our price: \$6,999.00



GE Dash 3000 Pro Patient Monitor with CO2
Compare price: \$9,999.00
Our price: \$4,999.00



Datex Ohmeda AS/3 Compact Patient Monitor
w/ 5 Agent Monitoring
Compare price: \$9,999.00
Our price: \$4,499.00



Criticare CSI 8100 Poet Plus Vitals Signs
Patient Monitor with Co2
Compare price: \$9,999.00
Our price: \$4,999.00



Datascope Passport 2 & Datascope Gas
Module II
Compare price: \$9,999.00
Our price: \$4,599.00



Invivo M12 3550A Anesthesia Monitor
Compare price: \$9,999.00
Our price: \$4,699.00



GE Datex Ohmeda Cardiocap/5 Patient
Monitor
Compare price: \$6,999.00
Our price: \$4,199.00



Datascope Passport 2 Patient Monitor w/ Co2
 Compare price: \$6,899.00
Our price: \$3,800.00



GE DASH 4000 Patient Monitor
 Compare price: \$5,000.00
Our price: \$3,299.00



GE DASH 3000 Patient Monitor with Printer
 Compare price: \$7,999.00
Our price: \$3,399.00



Solar 9500
 Compare price: \$6,700.00
Our price: \$3,300



CardioDynamics BioZ ICG Hemodynamic Monitor
 Compare price: \$6,999.00
Our price: \$3,499.00



Datex Ohmeda AS/3 Anesthesia Patient Monitor Co2
 Compare price: \$6,500.00
Our price: \$3,399.00



GE DASH 3000 Patient Monitor
 Compare price: \$6,999.00
Our price: \$3,499.00



Larsen & Toubro Star 50 Patient Monitor
 Compare price: \$5,999.00
Our price: \$3,299.00



Datex Ohmeda AS/3 Patient Monitor CO2
 Compare price: \$6,500.00
Our price: \$3,299.00



GE Marquette Eagle 4000 Vital Signs Patient Monitor with ECG NIBP SPO2 CO2
 Compare price: \$5,900.00
Our price: \$3,299.00



Datex Ohmeda AS/3 Compact Patient Monitor
 Compare price: \$6,999.00
Our price: \$2,999.00



DATEX-OHMEDA CAPNOMAC ULTIMA Anesthesia Monitor
 Compare price: \$5,999.00
Our price: \$2,999.00



PROTOCOL PROPAQ ENCORE 204 EL
 PATIENT MONITOR PROPAK
 Compare price: \$5,999.00
Our price: \$2,999.00



GE MARQUETTE SOLAR 8000 PATIENT ICU
 CCU MONITOR SYSTEM
 Compare price: \$7,999.00
Our price: \$2,999.00



Invivo Millennia Multiparameter Patient
 Monitor
 Compare price: \$4,999.00
Our price: Call For Price



Welch Allyn Atlas Vital Signs Patient Monitor
 Compare price: \$2,999.00
Our price: \$2,999.00



GE Dinamap Carescape v100 - BP SPO2
 Temp- NEW
 Compare price: \$3,342.00
Our price: Call for Price



DINAMAP PRO 1000 Vital Signs Patient
 Monitor
 Compare price: \$4,900.00
Our price: \$2,699.00



GE DASH 2000 Vital Signs Patient Monitor
 EKG ECG NIBP SPO2
 Compare price: \$3,999.00
Our price: \$2,699.00



GE Dinamap Carescape v100 - BP SPO2-
 NEW
 Compare price: \$3,180.00
Our price: Call for Price



SPACELABS 90369 PATIENT MONITOR
 NIBP SPO2 ECG NIBP
 Compare price: \$5,200.00
Our price: \$2,499.00



Protocol Propaq CS 242 patient monitor
 Compare price: \$4,000.00
Our price: \$2,499.00



MDE ESCORT II + PATIENT MONITOR NIBP
 ECG TEMP
 Compare price: \$3,999.00
Our price: \$2,499.00



Datex Ohmeda AS/3 Patient Monitor
 Compare price: \$6,500.00
Our price: \$2,499.00



Philips MP5T Patient Monitor
Compare price: \$4,999.00
Our price: \$2,499.00



Spacelabs Medical 90367 Patient Monitor
Compare price: \$3,999.00
Our price: \$2,499.00



Agilent Philips M3 M3046A Patient Monitor
Compare price: \$3,999.00
Our price: \$2,499.00



GE Marquette Eagle 4000 Vital Signs Patient Monitor with ECG NIBP SPO2
Compare price: \$5,900.00
Our price: \$2,199.00



Novamatrix Capnogard ET CO2 Patient Monitor
Compare price: \$2,099.00
Our price: \$2,099.00